



What if waterline maintenance could be simple?

Sterisil® is the solution

“

My patients notice a bad taste with iodine.

“

Why do 1 of every 3 tests fail?

“

Clean and reinstall the bottle every day - such a pain!

“

Wish there was something simpler and more reliable.



Simple Maintenance

Be more efficient by doing less with Sterisil®

Our system drastically reduces the number of steps needed for DUWL maintenance, **reducing staff burden** and **increasing productivity**.

240	Daily Maintenance*	0
	Clean/dry bottle, wipe straw, reinstallation	
12	Monthly dummy straw**	0
	Installation, shock, reinstallation	
4	Quarterly dummy straw***	0
	Installation, test, reinstallation	
2	Replace straw****	1
	Installation, shock	
258	TOTAL ANNUAL STEPS PER OPERATORY	1
With Iodine ¹		With Sterisil

Simplify your water management with Sterisil® and focus on what matters most – *your practice and your patients.*




*Using 1 liter per working day
 **Monthly shock requires dummy straw with iodine
 ***Quarterly testing requires dummy straw with iodine
 ****2 straws per year

Simple Efficacy

Be more confident and pass more tests with Sterisil®

Sterisil® powerful non-oxidizing **Silver Ion** technology is **20x more effective** (≤ 10 CFU vs 200 CFU) ensuring water safety and compliance.



Sterisil® Straw starts, and remains above the worry-free zone



With Sterisil®, say goodbye to these problems:

- Waterline tests that fail more than 30% of the time.²
- Patient complaints about the foul taste and odor of your dental water.
- Worrying about waterline contamination risk with unstable iodine treatments.

Sterisil® Silver Superior DUWL Maintenance

Infection control made easy

Our Promise: 365 days of worry-free water

SHOCK Simply

Integrated shock for
hassle-free water treatment.



**Citrisil™ Shock
Tablets**

TEST with Confidence

Pass/Fail results in 15 minutes.



**FASTCheck15™
In-Office Dental Water Test**

MAINTAIN with Ease

EPA registered, **50x** lower
than CDC guidelines.



**Sterisil® Straw
(Municipal or Distilled)**

Choose Simple. Choose Silver. Choose Sterisil.



**Solmetex
is here to
assist**



My Solutions Center

A simple online resource to manage
all your dental water needs



Subscription Service

Make replenishing the
products you need easy



General Inquiries

**877.207.1551
Solmetex.com**

Contact your Solmetex Water
Specialist or Preferred Dealer for
a free water consultation

References

1. Crosstex International. (2015, December 15). DentaPure® Advanced, Easy to Use and Cost Effective Clinical Water Treatment Technology, Instructions For Use. PDF. Retrieved from https://www.crosstex.com/sites/default/files/public/educational-resources/technical-specifications/dp_ifu_6.pdf.
2. Study cited in "Treating and Monitoring Dental Water" Compendium of Continuing Education in Dentistry, by Nancy Dewhirst, RDH, BS & John A. Molinari, PhD. Data collected from 22,196 consecutive waterline tests (R2A, heterotrophic plate count tests) conducted in 2017 by ProEdge Dental Water Labs. All product types provided by dental offices. The efficacy chart baseline of effectiveness shows that an iodine two-straw system is less effective than one Sterisil silver straw.